

NANO Beam position monitor and feedback controller

Robert Russell, Julian Adams, David Barber, & Nigel Boulding




FMB Oxford




FMB Oxford NANO BPM

In-vacuum detector head, image processing unit (B100) and control computer (rack mount or desk top)

FMB Oxford exclusively licences this technology from the University of Manchester




The NANO BPM operates by scattering a small fraction of the beam from a low z foil (KAPTON, Mica or Diamond) which is then analysed by an in-vacuum pinhole camera



The NANO BPM comes with a full suite of control routines.



The system operates on a Windows 7 server. A LabView executable interfaces to the B100 and serves EPICS process variables.

FMB Oxford also supply a full EPICS database and EDM GUI set for integration.





The NANO BPM comes with four independent 0-10V feedback outputs. Each output can be driven by x or y peak position, centroid position, IO or FWHM.

The NANO BPM can be driven in and out of the beam with a linear drive, the in-position repeatability is achieved with a ball and cone





NANO BPM image and profiles of the I22 beam at Diamond Light Source





Analysis of the I22 beam in the vertical

Panel 1 is the vertical position of the peak in microns, the beam size is approximately 300µm the RMS drift is approximately 0.5µm up with a peak to peak maximum of 1.5µm. Panel 2 is the 3D reconstruction of the beam profile over time



Comparison of the resolution of the NANO BPM with the performance of a quadrant ion chamber at ESRF



Feedback control of a double crystal monochromator. Left panel a 1keV scan without feedback and Right with feedback. The first points are the monochromator moving to the start position of the scan.

Specifications

Operating principal	Low Z scatterer, pinhole focussed CMOS imager
External accuracy	< 100nm (200µm beam, 10 ¹² photons/sec)
Integration time	0.02 – 5 seconds
Digitization	16 bit
Power input	12V DC

Outputs

Number	Four
Signal Type	Analogue voltage 0-10V
Mapping	Independently mappable

Ordering Information

FMB Oxford NANO beam position system

AHD2402

FMB Berlin

Tel: +49 (0)30 67 77 30-0
info@fmb-berlin.de

www.fmb-oxford.com
www.fmb-berlin.de

FMB Oxford

Tel +44 (0)1865 320300
sales@fmb-oxford.com